



# Enovation of Hydrophobic IOLs

**enova**<sup>®</sup>  
Preloaded Hydrophobic Monofocal IOL



## Enovation of Hydrophobic IOLs



100%  
Glistening-Free Material




Outstanding  
Biomechanical Properties



Clinically Proven  
Low PCO Rate



Pre-Conditioning Free




Preloaded



# Enovation of Hydrophobic IOLs

Hydrophobic IOLs have evolved significantly over time, and Enova® represents the pinnacle of this evolution. That's why we call it "Enovation".




## enova® Generation

Enova® is the only 100% glistening-free<sup>1</sup> hydrophobic acrylic IOL that requires neither pre-hydration nor storage in saline. Enova®, which has 7% water content, is dry-packed and has excellent optical and mechanical properties with a preloaded system.



## 2<sup>ND</sup> Generation Hydrophobic IOLs

Hybrid polymers, which include less than 5% water content, demonstrate a so-called control over water uptake and improved resistance to glistening formation, albeit at reduced levels. However, IOLs made from such materials either have poor mechanical properties or require both pre-hydration and storage in saline<sup>2</sup>.



## 1<sup>ST</sup> Generation Hydrophobic IOLs

Acrylic IOLs, with less than 1% water content<sup>3</sup>, develop various levels of glistening post-implantation due to uncontrolled water intake into the IOL polymers.

1- Glistening Analysis in Enova® Hydrophobic Acrylic Intraocular Lenses / In-vitro Study Evaluating the Tendency of Different Intraocular Lenses to Form Intraoptical Glistenings by the University of Utah

2- Bausch & Lomb. enVista Directions for Use

3- Comparative analysis of in vitro accelerated glistening formation in foldable hydrophobic intraocular lenses. International Ophthalmology Tandogan, T., Auffarth, G. U., Choi, C. Y., Son, H.-S., & Khoramnia, R. (2021).




# Enovation of 100% Glistening-Free IOL Material

The Enova<sup>®</sup> IOL Material is the first 100% Glistening-Free hydrophobic acrylic IOL that does not require pre-hydration and storage in saline solution!

The Enova<sup>®</sup> PGF3 IOL is dry-packed and boasts exceptional optical and mechanical properties.


## Glistening Formation in IOL

Water molecules bind to certain chemical groups through weak hydrogen bonds. Over time, more water molecules diffuse into the polymer network and bind preferably to other water molecules, which forms clusters referred to as "glistening."



## ENOVA<sup>®</sup> 100% Glistening-Free IOL

The unique composition of Enova<sup>®</sup> material allows the uniform hydration of specific sites, controlled water uptake, and resistance to glistening formation.



unique hydrophobic  
IOL material

in-house production

EST. 2009



# Enovation of 100% Glistening-Free IOL Material

## Conclusion by the University of Utah



*In vitro glistenings study: University of Utah*

**Study:** *In vitro study evaluating the hydrophobicity of different intraocular lenses*

*In vitro intraocular glistenings*

**STUDY REPORT**

International Ocular Research Center  
John A. Moran Eye Center  
University of Utah

**Specialty:** VSV Biotechnology



*In vitro glistenings study: University of Utah*

**Conclusions:** Enova® hydrophobic acrylic intraocular lenses exhibited no glistening formation after hydration and variation of the temperature. Tecnis intraocular lenses exhibited trace glistening formation, and AcrySof intraocular lenses exhibited mild glistening formation in these in vitro test conditions. The new Enova® intraocular lenses showed no surface haze and glistenings when compared with other commercially available hydrophobic acrylic IOLs as AcrySof and Tecnis intraocular lenses.

Liliana Werner, MD, PhD

Nick Mamalis, MD

**Week 1:** The findings at this time point were generally similar to those on Day 1. Overall, whenever glistening formation was observed in this study at week 1, it was mostly within the central 4.0 mm of the IOL optic.

- Enova® IOLs: No glistening, no haze.
- Acrysof IQ IOLs: Mild optic haze (giving the lens a slight yellowish/brownish discoloration under light microscopy) and mild glistening formation. Diameter of the glistenings: 10 to 20 microns.
- Tecnis IOLs: Moderate central optic haze (giving the central part of the optic a yellowish/brownish discoloration under light microscopy) and trace glistening formation. Diameter of glistenings: 25 microns.

In an attempt to quantify glistening formation within the lenses, the number of glistenings or microvacuoles (MV) that were well focused in the X200 light photomicrographs (area of 0.35 mm<sup>2</sup>) were counted, and the results were converted to MV/mm<sup>2</sup>.

IOL	MV/mm <sup>2</sup> Week 1
Enova®	0
AcrySof IQ	8.7
Tecnis	2.9

Table 1 : Number of Microvacuoles Converted to MV/mm<sup>2</sup>



*In vitro glistenings study: University of Utah*

**Conclusions:** Enova® hydrophobic acrylic intraocular lenses exhibited no glistening formation after hydration and variation of the temperature. Tecnis intraocular lenses exhibited trace glistening formation, and AcrySof intraocular lenses exhibited mild glistening formation in these in vitro test conditions. The new Enova® intraocular lenses showed no surface haze and glistenings when compared with other commercially available hydrophobic acrylic IOLs as AcrySof and Tecnis intraocular lenses.


Liliana Werner, MD, PhD

Nick Mamalis, MD




## Outstanding Biomechanical Properties

Gentle and controlled unfolding process in the posterior chamber and no pre-warming or special pre-conditioning is required.



Introducing our groundbreaking IOL, delivering easy unfolding, special haptic design for great stability, and smooth injection capability.




Experience the convenience of controlled deployment for precise positioning and a seamless implantation process.



## Clinically Proven Low PCO Rate

Posterior capsule opacification (PCO) after cataract surgery is impacted by the intraocular lens' (IOL) design and material. Enova®'s new 100% Glistening-Free material minimizes the risk of PCO and Nd: YAG procedures after implantation.

In the multicenter studies performed on Enova® IOLs, PCO was evaluated on 320 eyes. After 1 year, the post-operative results showed that only 5% of the total eyes and 1% of total implantations had PCO, necessitating Nd-YAG laser treatment.



1- VSY Biotechnology Data on File, 2023.

2- RCOphth National Ophthalmology Database Audit Feasibility Study of Post-cataract Posterior Capsule Opacification 2021

## No Pre-Conditioning Required

A polymer's Glass Transition Temperature (Tg) is reached when the polymer changes from a rigid material to a soft material. Having a Tg of  $-2.0^{\circ}\text{C}$ , all IOLs with the unique Enova® material undergo a gentle and controlled unfolding process below standard operating room temperatures.

Thus, no warming or special pre-conditioning is required.

IOL	Tg (°C)	Glistening	Packaging State
enova®	-2.0	No	Dry
AcrySof Vivity®	15	Yes	Dry
Tecnis®	14	Yes	Dry





## Ready-to-Go Preloaded System

Due to its patented Rotaryjet technology, the Enova® Hydrophobic IOL with the Rotaryjet Preloaded IOL system provides a safe, efficient, and user-friendly delivery procedure. Its smart design allows for reliable surgery with smooth IOL implantation and reduced post-op risks.


### Step 1

Apply BSS,  
followed by the OVD



### Step 2

Rotate the mechanism  
through 90 degrees  
until you hear a  
'click'




### Step 3

Push the injector until the  
blue cushion is visible  
through the cartridge.  
Then, pull it back gently  
until it is automatically  
stopped by the barrier



### Step 4

Now it is ready to inject



## Technical Features

Enova <sup>®</sup> PGF3		
<b>Material</b>	Single Piece, 100% Glistening-Free, Hydrophobic Acrylic, Dry-Packed	
<b>Optic Design</b>	Monofocal, Biconvex Aspheric	
<b>Refractive Index</b>	1.53 (546 nm)	
<b>Glass Transition Temperature (Tg)</b>	-2°C	
<b>Water Content</b>	7%	
<b>Optic Diameter</b>	6.00 mm	
<b>Overall Diameter</b>	13.00 mm	
<b>Haptic Design</b>	C-Loop	
<b>Haptic Angle</b>	0°	
<b>Spherical Power Range</b>	Preloaded	From +6.00 D to +30.00 D (0.50 D increments)
	Standalone*	From 0.00 D to +32.00 D (0.50 D increments)
<b>Lens Color</b>	Clear	
<b>Photo Protection</b>	UV Filtration	
<b>Recommended Constants</b>	Ac A constant: 118.0 SRK-II : 119.03 SRK-T: 118.7 Haigis a0, a1, a2: 1.11, 0.4, 0.1 HofferQ pACD: 5.33 Holladay sf:1.55 Barrett Universall LF:1.73	
<b>Recommended Injector System</b>	Rotaryjet Preloaded System	Acrijetfly 2.2

\*Recommended injector system: Acrijetfly




enova®

Preloaded Hydrophobic Monofocal IOL

E.MNF.BRO.12/24.ENG\_rev04



VSY Biotechnology GmbH  
Esslinger Str.7 70771 Leinfelden-Echterdingen Germany  
contact@vsybiotechnology.com / www.vsybiotechnology.com

 vsy-biotechnology