

# Ocuva


ASPHERIC FOLDABLE INTRAOCULAR LENS



# Aspheric

**Ocuva**  
ASPHERIC FOLDABLE INTRAOCULAR LENS

- Aspheric
- Single Piece
- Foldable
- Hydrophilic Acrylic
- Square Edge
- UV Blocking



Perfect Image Quality  
With Ocuva Aspheric



|                | Ocuva A625 | Ocuva AB625         |
|----------------|------------|---------------------|
| Optic Type     |            | Aspheric            |
| Optic Size     |            | 6.00 mm             |
| Optic Design   |            | Biconvex            |
| Optic Material |            | UV blocking acrylic |
| Water Content  |            | 26%                 |
| Haptic Size    |            | 12.50 mm            |
| Haptic Angle   |            | 0°                  |
| Haptic Design  | Modified C | Balance Design      |

| A Constant         |       |
|--------------------|-------|
| Rec.Ac. A Constant | 118.0 |
| Rec.Op. A Constant | 118.5 |

| Diopter     |   |
|-------------|---|
| Power Range | From 0.00D to +32.00D (with 0.50D increments) |

| Refractive Index 546 nm |            |                       |
|-------------------------|------------|-----------------------|
| Dry                     | 20°C /35°C | 1.510 / 1.510 ± 0.002 |
| Wet                     | 20°C /35°C | 1.463 / 1.463 ± 0.002 |

| Refractive Index 589 nm |            |                       |
|-------------------------|------------|-----------------------|
| Dry                     | 20°C /35°C | 1.508 / 1.508 ± 0.002 |
| Wet                     | 20°C /35°C | 1.461 / 1.461 ± 0.002 |

| Barcode |   |
|---------|---|
| EAN     | EAN 13  |
| 2D      | All Ocuva series intraocular lenses are tracked by 2D Barcode system showing all the production process |

**AcriJET** *green* Injector and Cartridge  
CARTRIDGE & INJECTOR SYSTEM

Ocuva can be implanted with Acrijet Green 2.0 and Acrijet Green 2.2 Injector and Cartridge system.

| Injector          | Incision Size | Optic Type   | Maximum Diopter Range |
|-------------------|---------------|--------------|-----------------------|
| Acrijet Green 2.0 | ≥2.0 mm       | Ocuva A 625  | sph 25.0 D            |
|                   |               | Ocuva AB 625 |                       |
| Acrijet Green 2.2 | ≥2.2 mm       | Ocuva A 625  | sph 30.0 D            |
|                   |               | Ocuva AB 625 |                       |

